



B E S T O F T H E B E S T

www.becasmedikal.com.tr



We love to produce technology for the world!

www.becasmedikal.com.tr



We added our technology and experience on all of our devices as we took this road. We have made our studies and researches perseveringly in order to take better and firm steps forward.

Becas Family concept became our milestone. Maybe we didn't sleep for the success of our family members who believed in us for nights. We proceeded confidently. Because we knew that the Becas family deserves the highest quality of technology.

I offer my endless thanks to every point that believes in us. The world must know the Becas family, it is time to represent our country.

Yours truly.



Burcu Soydan
Becas Medikal CEO

HISTORY

Becas Medikal entered the sector in 2006 and started IPL production in 2008. It launched in 2010 the sales of the Dr. Rimpar products in Turkey market. Becas has achieved great success and added Isabelle Lancray's products in 2011 in its portfolio.

In 2013, it began with the sales of NPM International. In 2014, the company made a big move and started to produce diode lasers. In 2015, it started also to produce permanent makeup devices.

In 2016, with the understanding that attaches importance to congresses and seminars, it held the biggest Greetings Congress in Europe. In 2017, the company got the first place in the export in its sector. In 2018 B8 Diode Laser Device became the most preferred device of Turkey.



Certificates



Becas Medical, proving its successful business partnerships and sales success with its network in 14 different countries around the world, continues to grow by developing its technology and product range.

EUROPE

ASIA

AFRICA

Morocco
Algeria



The B8 Platinum is now waiting for you so that you can make applications
in four seasons and get excellent results!



TOUCH SCREEN



35 MINUTES PROCESSING TIME



- 10 ° C ICE CAP



4 SEASONS USE



10 SHOTS IN 1 SECOND




STERILIZATION FEATURE



EMERGENCY BUTTON





Before Application


After Application

> American Bar Cap

The American Bar Cap used in the B8 Platinum allows up to 40,000,000 shots due to product quality. It works much more efficiently and stably than the caps used in other devices.


> Quartz Sapphire Glass

Quartz Sapphire Glass prevents the hairs from sticking to the cap during the application and increases the comfort for the user and the customer.



> Ice Cap

The ice cap, which can cool down to -10°C , reduces the temperature to which the skin is exposed to normal level, minimizing the risk of burning.




> 10 Shots in 1 Second

The B8 Platinum is able to shoot up to 10 shots in 1 second with its adjustable program.



> Custom Software Support

Perfect results in the most successful and short sessions thanks to software for men and women!




> Cap Sterilization System

B8 Platinum offers a standard sterilization system that gives importance to human health. This patented and unique system is offered as a standard feature to its user in B8 Platinum.



> New Generation Cooling System

B8 Platinum, which succeeds to exceed the cooling standards, maximizes the overall cooling performance of the device, providing long-term and safe use.



> Cap Cleaning Slot

Offered as the first and only in the world, this patented feature is self-cleaning with the release of the B8 Platinum cap.

> Emergency Button

Stop the entire system with one touch in any emergency case. A system that allows you to use!

> 10.4" Touch Screen

The B8 Platinum features an 10.4-inch upright touchscreen integrated in the housing for quick and convenient use.

> Carbon Coating

The carbon coating used in Porsche vehicles makes the B8 Platinum case look stylish.



B-Eight Platinum technology blended with experience, the latest point of a diode laser system B-Eight Platinum is now with you.

Its double cooling system provides the coldest cap on the market with -10 C ice cap. This provides comfortable and painless application

Thanks to the special software for women and men in the B-Eight Platinum devices, it is possible to achieve the best results in the least number of sessions.

The device can sterilize its own cap so that it can be safely operated in different regions and different customers.

It is a privilege to work with the B-Eight Platinum ironing system, which can provide up to 10 shots per second, allowing a complete body treatment time of 35 minutes.



Usable in 4 seasons

Beight DIODE LASER

A New Miracle in Every Shot, Be Witness This Miracle!



EMERGENCY BUTTON



10 SHOTS IN 1 SECOND



35 MINUTES PROCESSING TIME



STERILIZATION FEATURE



- 5 ° C ICE CAP



TOUCH SCREEN



10 shots in 1 second

The B-Eight Laser Hair Removal Device allows the user to shoot up to 10 shots per second.



- 5°C Ice Cap

With the refrigerant in the cap, it goes down to -5 ° C and this shows the importance given To human health by reducing the burning risk to a minimum level.

Trigger System

The trigger system, which is very convenient for the user, provides an easy and comfortable application. The trigger button, which makes it easy to control the B-Eight, is always at hand!




- The B-eight diode laser is patented in the world with its patented sterilization feature that ensures the removal of all microorganisms from all living forms.
- Beight, which has 3 ready-to-use programs, performs and stops as many as the number of shots in the selected single, double or triple shots.
- In 2017 we ranked to the 1 place in export with our international experience. Beight diode laser device represents our country in the world arena.

B-eight Diode ironing laser hair removal became the brand in our country in 2015 and the laser hair removal applications were moved to another dimension. With the demand, the B-eight device has become very popular in the past years, now with the original ironing hair removal technology there is no stopping in the summer months and four seasons laser epilation can be performed.

Unlike other devices, the contact cooling cools down to - 5 ° C and the cap reaches a consistent icing level during application.

That's where it's painless and ache less. Treatment session times now take an average of 35 minutes for the whole body because it is able to shoot 10 shots per second and give you 12x12 spot of application cap. Patented sterilization technology eliminates the risk of infection.

- Featuring an 8-inch upright touch screen, the B-Eight provides fast and reliable application.
- It provides a fast and reliable application with 808 Nanometer with American Diode cap, 10 Hz speed, Peltiel and Climatic cooling.



Beight DIODE LASER

Laser Type	DIODE LASER (Ironing)
Power	600 WATT
Wavelength	808 nm
Spot Area	12 x 12 mm
Screen Size	8 '' Touch Screen
Cap Cooling	Peltier Cooling
Device Cooling	Climatic Water Cooling
Operating Voltage	220V / 15A max 50/60 Hz
Energy	SHR 12J / cm2



Beight URBAN

Life is Now More Colorful More Free!

Beight Urban colors the Life Even More!



B-eight URBAN



EMERGENCY BUTTON



10 SHOTS IN 1 SECOND



35 MINUTES PROCESSING TIME



TOUCH SCREEN



- 3 ° C ICE CAP



4 SEASONS USE



Body Color Option

You can choose one of the different body color options according to your business structure or taste.

An unrivaled technology that delivers 10 shots per second and the cooler at the tip of the cap goes down to - 3 ° C. Ironing Hair removal is a privilege and the B-Eight Urban diode laser system is the number one brand of this privilege.

- Featuring an 8-inch upright touch screen, the B-Eight Urban provides fast and reliable application.
- It provides a fast and reliable application with 808 Nanometer with American Diode cap, 10 Hz speed, Peltier and Climatic cooling.



Feel the Change in Aspect!





Cavitation Cap

The cavitation cap is used to render the applied area of the body in a painless, leak-free manner.

Display System

Thanks to its double cap, it provides the opportunity to see all the information thanks to the display screens in the application and caps in 2 different regions at the same time.



Small Cap

Thanks to its small cap, it provides more comfortable application in double chin and arm areas.



360° Cap Design

The most important feature that distinguishes the B-Cryo 360 ° from other devices is that it provides a more comfortable and accurate use thanks to its 360 ° winding and freezing of the applied area.





Touch Screen System


10" touch screen and easy-to-understand interface provides ease of use to the user.



Changeable Cap System

The 360 ° feature of the B Cryo 360 ° is available as standard in all three easily changeable caps. The B Cryo 360 ° shows its difference from its rivals with its 3 changeable caps.





OTHER BRANDS

Cooling Only with
Right and Left Peltiers

B Cryo 360 ° Cap Structure

It is 360 ° cooling of the vacuumed area. In all devices except the B-Cryo 360, this cooling is not 360 °, but only two metal peltiers on the left and right inside the cap. That is why the B-Cryo 360° is the only device that can provide a balanced and maximum thinning in the treated area. B-Cryo's 360 ° cooling cap minimizes cold burns and all other risks and maximizes success.

Cap Structure of Other Brands

Other cooling caps used in the industry, only the right and left of the peltier contacted by making contact in the area prevents a healthy operation. Increasing the risk of cold burns to a high level and therefore the success rate is very low.



Always Worth Being Young

You Will Not Believe in the Rejuvenation Effect of the Golden Needle

> Vacuumed Special Cap (Patented)

Thanks to the vacuum fractional radio frequency technology, radio current is delivered to the target points under the skin with gold-plated needles without any pain or sensation. Thus, the skin self-repair and renewal mechanism is activated. There are three different needle options with 10, 25, and 64, the expert can choose according to the size of the application area.


> Cooler Cap

Immediately after the fractional radiofrequency application, the cryo cap is applied, thus preventing the feeling of redness and pain completely and providing a comfortable application.

> Smart Display

As it is the case in almost all our devices, it brings the user with the technology together with its easy interface and touch screen in order to perform their sessions easily for the user.





> Different Needle Options


With the replaceable cap system of B-Gold, gold needles can be used in 3 different sizes as 10 pieces, 25 pieces and 64 pieces. There is also a needlefree radio frequency tip.

> Vacuum Cap System

Special patented vacuum cap provides comfortable application without pain and ache.

> Gold Plated Needles

B-Gold's gold-plated needles do not harm the skin in any way with different depth adjustment






Usage Areas

- » Eliminating fine wrinkles,
- » Recovery of sagging in the body,
- » Skin rejuvenation,
- » Acne scarring, pitting and scar removal,
- » Treatment of burn scars,
- » Treatment of bruises
- » Cracks in pregnancy,
- » Pore tightening,
- » Acne treatment,
- » It plays an important role in recovery of neck and décolleté.

Gold Needle, Radio Frequency, specially
produced for Dermapen applications





For Skins Like Pure Water

Hydra, with 4 different solutions used with caps deep cleansing is applied to each skin type and the peeling and blackhead cleaning can be done without touching. The device is distinguished from the devices in the sector with its self-cleaning feature.

Deep cleaning, disinfection and tightening of the pores are ensured with the frimotor and high frequency heads, while the SOS-Serum oxygen spray head ensures that the skin is saturated with oxygen, the cell renewal is activated at the maximum level and the product penetration to the deepest depths.

Led therapy provides collagen and elastin stimulation. And with the Gemeo capsule application technique only available in this device, it provides instant visible results in skin spots and wrinkles. Oxygenation of the skin is increased. It is a very privileged and personalized application. Thanks to the meomassage feature, it enables the acceleration of micro circulation.



Camera and HD Screen

Thanks to the special camera and HD Screen of the device, the skin is analyzed in detail and healthy and the correct care program is determined before starting the care. At the same time, the camera is very useful in performing hair root analysis before laser epilation. It was observed that the clients who saw the image of their skin on the screen with 500 times enlargement were observed to be much more motivated and disciplined in skin care and product use.

Technical Specifications

Dimensions	48 x 40 x 110 cm
Color	White
Net Weight	24 Kg.
Output Voltage	50 VA
Equipment Safety Class	I Class BF Type
Solution	4 Different Options
Cap Option	6 Different Options

Aquafacial Handle

With this special cap user with 4 different solutions, lactic acid, salicylic acid or hyaluronic acid can be given to the skin according to the need, and its withdrawn with vacuum to deeply cleanse and moisturize the skin.

Different size tips are available for each region. Furthermore, the device can clean itself with the 4. solution. This is a very important feature not found in other Hydrafacial devices.

Oxygimeo Handle

Depending on the skin needs of the person, there are two types of tips that can be selected related to stain or antiaging. This tip is attached to the Oxygimeo cap of the device and after applying the appropriate gel to the skin, it is performed on the surface of the skin with this cap.

Skin cell metabolism, which is intensely exposed to carbondioxide, is saturated with oxygen more than usual, spots are removed and wrinkles disappear.

Frimator Brush Head

Thanks to this head, which can be used with 5 different tips, the dead cell layer on the skin is deeply cleaned and the skin breathes. With this application, which is indispensable for a perfect skin care, dry and / or wet peeling process is performed and blood circulation is accelerated. The rotation speed of the brush can be adjusted according to the condition of the skin.

High Frequency

Used with probes in 4 different sizes and shapes, this head allows better absorption of the active substances applied by the skin, while at the same time helping to narrow the pores and get rid of acne. It is germicidal and antibacterial due to the ozone formation. It opens clogged pores, allows the skin to breathe and increases blood circulation, allowing the skin to gain a healthier appearance.

Serum Oxygen Spray Gun

This head, which allows dermocosmetic products to be penetrated under the skin with high pressure, is an alternative to needle methods. It provides the skin to saturate and moisturize with oxygen. It revitalizes the skin, accelerates blood circulation, and is an indispensable part of antiaging and spot treatments.

B-RF RADIOFREQUENCY


For being rejuvenated stay on that frequency!
Youth Effect on Face with B-RF Radiofrequency


Unlike conventional radiofrequencies on the market, the B-Rf Radiofrequency Vacuum provides three different sizes of caps that allow the radio waves to reach the subcutaneous target points more easily.

Radiofrequency technology enables the localized fats that are resistant to diet and exercise to be disintegrated by heat and discharged through the lymph channels. Vacuum nozzles increase the blood circulation in the application area thanks to 4 different alternatives which can make vacuum continuously and occasionally, the elastin and muscle tissue are stimulated to achieve tremendous results in both regional thinning and skin rejuvenation.


B-Rf radiofrequency treatments;

- Antiaging
- Spot Treatment
- Regional Thinning
- Cellulite
- Collagen and Elastin Bond Repair







VACUUM TECHNOLOGY




3 CAP OPTION



MULTI POLAR FEATURE



4 SEASONS USE



EASY CLEANING



Technical Specifications

Power	400w
Radio Frequency	0.8 MHz (Multi-polar RF)
Vacuum Range	318mmHg
Dimensions	(WxDxH mm) 320 x 451 x 330
Weight	8.5 Kg
Electrical requirement	220/260 50-60hz

SKIN ANALYSIS DEVICE



As a result of this analysis, the beautician may have a healthy idea about the right skin care alternative that the customer needs. The customer who sees this analysis on the screen can better understand the importance of skin care as she/he can see the details under the skin that she/he cannot normally see with the naked eye. According to the analysis, the device also recommends DR. RIMPLER products that the customer should use according to the skin type.

Thanks to our device, which provides the opportunity to compare the analyzes performed on different dates, both the customer and our beautician can follow in detail the positive effect of the treatments on the skin. The analyzes can be recorded in the device memory. It is one of the essential devices in terms of bringing new skin care customers to salons and supporting product sales.

PERMANENT MAKE-UP DEVICE

- » Touch screen
- » Multi-Processing with One Device
- » 17 Kinds of Needle Scales
- » Patented Roller System in Hair Simulation
- » Medical Pigmentation
- » Permanent Makeup and Tattoo Removal
- » Dermapen
- » Automatic Amperage Adjustment Special to the Processing

Medical Mode —○

Eyebrow —○

Eye —○


Lip —○





Resolving Skin Problems Now in Your Hands !

- Beauty Spots and Wart Treatment
- Scar Treatment
- Eyelid aesthetic
- Cracks Treatment
- Spot Treatment
- Removing Lip Upper Lines
- Abdominal Tightening
- Wrinkle Removal
- Face and Neck Lifting
- Acne Treatment
- Post-Acne Scar Treatment



The face on mole



Face spotted



Warts



CALL CENTER

+90 212 885 99 99

www.becasmedikal.com.tr | info@becasmedikal.com.tr

PRODUCTION

Becas Medikal A.S.

Halkalı Merkez Mahallesi Dereboyu Caddesi Balcon Plaza Kucukcekmece / Istanbul

